


# PEFC : Aberration corrected STEM-EDX analysis of core-shell nanoparticles.

Clear elemental maps of small nanoparticles can be acquired by combination of aberration correction, high sensitive EDX analysis, and low accelerating voltage measurement (80kV).  
 Distribution of Pt shell thickness in Pd-Pt core-shell structure can be analyzed precisely.


## Aberration corrected STEM-EDX analysis



Clear elemental maps of small nanoparticles about 5nm in size can be acquired by combination of aberration correction, high sensitive EDX analysis, and low accelerating voltage measurement (80kV).

## Precise analysis of core-shell structure


- Analysis procedure —
- ① Preparing profiles of the number of Pd and Pt atoms with different shell thickness.
  - ② Estimation of effective STEM probe size and convolution it to the profiles.
  - ③ Comparing with experimental results.



The number of shell layers is determined in an atomic level by fitting core-shell models and experimental EDX profiles. \*

\* S. Inamoto, Y. Otsuka, K. Kobayashi (Daido univ.), M. Hori (Daido univ.),  
 The 69<sup>th</sup> Annual Meeting of the Japanese Society of Microscopy.

## Evaluation by STEM image simulations



Results of the HAADF-STEM image simulations show the 3-layer model is similar to the experimental data.

Sample provided by Brookhaven National Laboratory