


# Online Accurate mass measurements of evolved gas ~TPD-TOFMS analysis~

Toray Research Center, Inc. established the TPD-TOFMS\* apparatus, which consists of a temperature-controlled furnace with a time of flight high-resolution mass spectrometer. TPD-TOFMS is applied to determine evolved gas species. Based on the advantage of TOFMS, it is possible to distinguish components of same integer mass number (e.g. CO or N<sub>2</sub>, SO<sub>2</sub> or S<sub>2</sub>).

## Equipment overview



Each evolved gas species can be detected as a function of heating time or temperature.

\* Temperature Programmed Desorption-Time Of Flight Mass Spectrometry

## Equipment specifications


Temperature Range : R.T.~1000°C  
 Heating rate : 1~50°C/min  
 Atmosphere : He, O<sub>2</sub>/He etc.  
 Mass range :  $m/z$  2-300  
 Mass accuracy :  $\pm 0.001$

※ Depending on experimental Condition

## Measurement of a material for sulfide all-solid-state batteries\*※

\* Sample was provided by Dr. M.Tabuchi and Dr. T.Kojima in AIST, Japan.


Outgas safety assessment and degradation in heating process are investigated via evolved gas analysis.




## Experimental Condition

Temperature Range : R.T.→500°C, 10°C/min  
 Atmosphere : He  
 Mass range :  $m/z$  10-300  
 ※ The sample was transported to apparatus without exposure to air.

☞ Monitorable evolved gas analysis as well as conventional TPD-MS. Moreover, it is possible to distinguish outgas components of same integer mass number via high resolution of mass detection.



☞ By analyzing the mass spectral data, S<sub>2</sub> peak (Fragment ion of cyclic S<sub>6</sub>, S<sub>8</sub>) are identified near SO<sub>2</sub> peak.



☆ Based on the TOFMS detector, we can separate SO<sub>2</sub> and S<sub>2</sub>. The precise interpretation of thermal degradation (e.g. progress of oxidation and/or change in chemical structure) can be discussed by TPD-TOFMS analysis.