


# Deep Proteomics

Here we demonstrate a new LC-MS/MS analysis method for peptide mixtures with high resolution by combining ion exchange chromatography and C18 long monolith column chromatography. Analysis of standard cultured cells resulted that more than 6000 proteins were identified at a quantifiable level.


## Deep proteomics workflow

All biological samples can be analyzed




## Analysis for human cultured cell (HeLa cell)

### Effect of SCX and monolith column




The number of identified proteins was significantly increased.

### Reproducibility of quantitation



Correlation factor : >0.999

### Mapping to signaling pathway



Many identified proteins were mapped to the signal pathway.

### Examples of the analyzed data

試料1	正規化したピーク強度						MassPro	UniProt	蛋白質名		
	測定1	測定2	測定3	測定1	測定2	測定3					
1.95E+08	1.35E+08	1.27E+08	1.66E+08	1.47E+08	1.38E+08	1.89E+08	1.71E+08	1.22E+08	301	COX19B	Mitochondrial-associated protein RPE65 family member 3
1.98E+08	1.35E+08	1.27E+08	1.66E+08	1.47E+08	1.38E+08	1.89E+08	1.71E+08	1.22E+08	302	COX19B	Mitochondrial-associated protein RPE65 family member 3
1.42E+08	9.58E+07	8.61E+07	7.73E+07	6.56E+07	7.19E+07	1.37E+08	1.15E+08	1.10E+08	303	PI1137	Mitochondrial-associated protein 2
2.82E+08	1.98E+08	1.20E+08	1.55E+08	1.80E+08	1.38E+08	2.81E+08	1.41E+08	1.73E+08	131	COX19C	Protein NDRC2
5.98E+08	4.13E+08	3.95E+08	5.49E+08	5.28E+08	4.81E+08	5.70E+08	4.78E+08	4.66E+08	65	PI1988	Mitochondrial-associated growth factor
1.91E+08	1.01E+08	6.83E+07	1.88E+08	1.23E+08	1.23E+08	3.30E+08	1.98E+08	1.37E+08	62	PA6842	Dyogenin synthase kinase 3 alpha
6.98E+08	6.03E+08	3.92E+08	8.43E+08	4.38E+08	3.23E+08	1.33E+09	5.21E+08	6.94E+08	424	PA6841	Mitochondrial-associated protein 1b
4.65E+08	3.47E+08	1.82E+08	5.73E+08	3.92E+08	3.03E+08	7.66E+08	3.96E+08	3.88E+08	701	PO4792	Heat shock protein beta-1
2.21E+08	2.94E+08	3.80E+08	1.92E+08	1.83E+08	2.23E+08	4.16E+08	3.23E+08	3.98E+08	313	PI10443	Biphenyl
6.27E+08	5.91E+08	5.77E+08	5.32E+08	4.66E+08	4.51E+08	8.28E+08	6.68E+08	6.98E+08	432	2SRI410	Biphenyl-adenosine
8.28E+08	4.98E+08	3.91E+08	7.67E+08	3.77E+08	3.52E+08	9.04E+08	5.98E+08	3.72E+08	310	PI1937	Insulin-like growth factor 1
7.97E+08	5.91E+08	3.75E+08	6.71E+08	5.18E+08	3.71E+08	6.88E+08	5.39E+08	4.41E+08	583	PI1137	Mitochondrial-associated protein 2
5.17E+08	3.59E+08	3.13E+08	5.86E+08	5.18E+08	3.00E+08	6.98E+08	4.22E+08	4.20E+08	944	COX14P0	Protein NDRC4
5.26E+08	4.90E+08	3.78E+08	5.12E+08	3.82E+08	3.86E+08	7.31E+08	4.95E+08	4.93E+08	351	COX19D2	Mitochondrial-associated protein 12
3.68E+08	3.11E+08	2.16E+08	3.68E+08	2.42E+08	2.27E+08	4.81E+08	4.07E+08	3.42E+08	794	PI19559	Mitochondrial-associated protein 1A
6.94E+08	2.94E+08	1.16E+08	5.26E+08	2.27E+08	1.34E+08	6.91E+08	2.14E+08	1.75E+08	144	2SRI403	Protein NDRC2
2.08E+08	1.73E+08	1.36E+08	1.63E+08	1.50E+08	1.30E+08	2.18E+08	1.74E+08	1.64E+08	561	2SRI408	Mitochondrial-associated protein tau
4.88E+08	4.43E+08	3.79E+08	4.77E+08	4.08E+08	4.57E+08	5.27E+08	5.32E+08	3.97E+08	310	COX19B	Dihydropyrimidinase-related protein 2
4.08E+08	3.18E+08	2.92E+08	3.45E+08	2.86E+08	2.74E+08	4.60E+08	3.51E+08	3.06E+08	110	CI1075	Beta-thrombomodulin protein (MIM: 61241)
4.48E+08	3.43E+08	2.88E+08	5.09E+08	3.86E+08	3.03E+08	7.59E+08	3.73E+08	3.30E+08	671	COX19B	Mitochondrial-associated protein RPE65 family member 3
6.88E+08	3.81E+08	1.94E+08	4.47E+08	2.86E+08	2.16E+08	6.90E+08	3.26E+08	2.18E+08	584	2SRI478	Biphenyl
8.84E+08	5.83E+08	3.02E+08	5.99E+08	5.11E+08	4.12E+08	9.22E+08	7.43E+08	6.11E+08	367	2SRI410	Biphenyl-adenosine
3.21E+08	2.46E+08	1.73E+08	3.09E+08	2.86E+08	2.16E+08	3.32E+08	2.17E+08	2.16E+08	313	PI1938	Biphenyl-derived growth factor
7.54E+08	5.40E+08	3.08E+08	6.87E+08	6.25E+08	4.33E+08	7.07E+08	5.18E+08	4.01E+08	317	2SRI478	Dihydropyrimidinase-related protein 2
2.55E+08	1.81E+08	1.29E+08	1.77E+08	1.25E+08	1.07E+08	2.44E+08	1.46E+08	1.44E+08	714	2SRI403	Protein NDRC2
2.14E+08	1.58E+08	1.36E+08	1.82E+08	1.28E+08	1.12E+08	2.57E+08	1.51E+08	1.57E+08	712	2SRI478	Biphenyl cluster assembly enzyme (MIM: 61240)
1.14E+08	7.80E+07	6.71E+07	1.11E+08	9.92E+07	6.07E+07	1.83E+08	1.26E+08	1.19E+08	407	COX19E7	B-ubiquitin-protein ligase C14P
1.25E+08	6.94E+07	1.15E+08	5.91E+07	3.26E+07	1.17E+08	5.30E+08	2.48E+08	1.92E+08	691	2SRI408	Mitochondrial-associated protein 1
6.91E+08	6.79E+08	1.99E+08	8.66E+08	4.26E+08	3.40E+08	6.69E+08	3.91E+08	3.17E+08	414	2SRI408	Mitochondrial-associated protein 1

We can provide data suitable for statistical analysis. We also support pathway analysis for the data.