


# Highly sensitive analysis of odorous component in gas

The online concentration system combined with GC can transfer odorous volatiles (number of C in Alkane > 2) desorbed from a sample into the GC column, and consequently can lower detection limit. The application to estimating odorous components desorbed from sulfuric compounds is presented.

## Analytical tech. of gas



Sampling method	Analytical method	Target components
Sampling bag (non concentration)	GC, GC/MS, FT-IR	H <sub>2</sub> , N <sub>2</sub> , O <sub>2</sub> , CO, CO <sub>2</sub> , NO <sub>x</sub> , SO <sub>x</sub> , Organic components (C1~8)
<b>Online concentration -GC/MS</b>		Organic components (C 2~13)
Adsorbent (Offline concentration)	GC/MS	Organic components (C4~30)
Solution (concentration)	GC/MS, HPLC, Ion Chromatography	Organic components (C10~) Acid or basic gas (NO <sub>x</sub> , SO <sub>x</sub> , etc.) Aldehyde, Ketone(Derivatization)

Detection sensitivity is 20 times higher (with n-Butane)

## Application

Sulfuric compounds

### Online concentration-GC/MS analysis

of gases desorbed from sulfuric compounds

Condition of desorption • Heating temp. : R.T ~ 200°C

Atmosphere : Air, N<sub>2</sub>, etc.

Table Qualitative result of peak No.1~4

Peak No.	Components	Odors*
1	Sulfur Dioxide (SO <sub>2</sub> )	The irritating odor such as the egg which became added.
2	Aldehyde	The irritating odor that was burnt and sweet-sour
3	Thiol	The skunk's intense acrid odor
4	Amine	The odor such as ammonia

\*Reference : SDS, Handbook of Offensive Odor Control Law (in Japanese), etc.




Fig. GC/MS-Total ion chromatogram

SO<sub>2</sub>, aldehyde, thiol, amine, etc. are estimated as the source of odors from Sample B.

## Other applications

Fluoric compounds

Nitrile, aldehyde, etc. which contained fluorine are estimated as the source of odors from the sample

Plastic products

Unsaturated hydrocarbon (C5), aldehyde, thiazole, etc. are estimated as the source of odors from the sample



Toray Research Center, Inc.