
## GC/MS Analysis under GMP — Highly Sensitive Analysis for Impurities in Nucleic Acids

In the quality test of nucleic acid medicines, plenty of samples can not be used for analysis, unlike the case of small molecule drugs. Especially in the analysis for impurities, highly sensitive analysis is need to perform the test using limited amount of valuable samples. In this report, the residual solvents contained in a commercially available nucleic acid sample were measured by highly sensitive headspace–GC/MS.


## 1. Quality Test Items for Nucleic Acid Medicines





SOP: Operation and maintenance of headspace-GC/MS (SOP No. NME18500)

## 2. Measurement of Residual Solvents in the Nucleic Acid by Headspace-GC/MS




The vapored sample is absorbed and concentrated at the trap tube, and introduced to GC/MS.

The sample sealed in a vial is heated in an oven and vapored. The vapor is introduced to a trap tube or a loop.

<Figure of a Headspace Device>

Improvement of detection capability for dichloromethane (a typical solvent used for nucleic acid synthesis)




Residual dichloromethane in artificially synthesized nucleic acid was measured by the trap mode.

Sample: Artificially synthesized nucleic acid

(C-Myc antisense DNA, sulfurized)

Sequence: CACGTTGAGGGGCAT Sample concentration: 1 mg/mL



Chromatogram of the artificially synthesized nucleic acid by the trap mode

A trace amount of dichloromethane (4ppm for the synthesized nucleic acid) was detected.

- S/N ratio of dichloromethane increased by 30 times in measurement by trap mode headspace-GC/MS.
- Quantitation of dichloromethane was performed using very small amount (1 mg) of the nucleic acid.

