

# Distribution of Impurities in semi-conductor device by NanoSIMS

NanoSIMS is specialized technique for imaging with high lateral resolution of  $\sim 50\text{--}200\text{nm}$  with high sensitivity. Furthermore, depth profiling can also be conducted using primary beam with higher current density enough to sputter up to an interested depth of sample.

## NanoSIMS 50L


### <Instrument>



(AMETEK HP)

Multi-collection system

### <Configuration>




### <Specification>

Primary ion species	$Cs^+$ , $O^-$
Beam size	$O^-$ , $Cs^+$ : $< 50\text{ nm}$
Detection limit	ppb $\sim$ ppm
Type of mass analyzer	Magnetic sector type
Number of detected elements	7
Analytical depth	10 – 500 nm

- ◆ High lateral resolution of 50 nm
- ◆ High sensitivity
- ◆ High mass resolution


## Imaging of Al and P of SiC-MOSFET

Al




Al includes residual component at the surface.

P



## Depth profiles of implanted elements in Si substrate



### <Measurement condition>

- ✓  $Cs^+$  beam, 16 keV
- ✓ Beam size :  $> 1\ \mu\text{m}\phi$
- ✓ Raster area :  $10\ \mu\text{m}\square$
- ✓ Detected area :  $1\ \mu\text{m}\square$

### Detection limit or Background level (atoms/cm<sup>3</sup>)

**C** :  $1.5 \times 10^{18}$   
**P** :  $1.5 \times 10^{16}$   
**As** :  $4.0 \times 10^{16}$